


Disposable CPAP filter for Philips DreamStation 2

Item No.: HC-139



Product Function:

CPAP filters are designed to reduce airborne contaminants such as dust, pollen, and microorganisms from the intake air, ensuring that clean and filtered air is delivered to the user during CPAP therapy.

Compatibility	Philips DreamStation 2
Classification	Class I according to (EU) 2017/745, Annex VIII
Dimensions & Tolerances	45.0 × 21.0 mm
Color	Filter medium: white Frame: gray
Filter Type	Disposable
Outer Material	Polypropylene
Filter Media	Electrostatic non woven cloth
Packaging	12 pcs per pack
Filter Efficiency	99.9% (NaCl, 0.3 μm @ 5.33 cm/s, ±2%)
Airflow Resistance / Pressure Drop	Tested according to ISO 29463
Performance Criteria	Maintains structural integrity under humid conditions. Does not degrade or shed particles during use.
Test & Validation	Pressure drop value 4.2 pa
Sterility	Non-sterile
Storage & Transport Conditions	Temperature: +5 °C to +38 °C Relative humidity: < 75%, non-condensing
Shelf Life	3 years from the date of manufacture when stored in a sealed bag and protected from direct sunlight