


Hypoallergenic CPAP-Filter for Resmed Airsense™ S9 & 10

Item No.: HC-151



Product Function:

CPAP filters are designed to reduce airborne contaminants such as dust, pollen, and microorganisms from the intake air, ensuring that clean and filtered air is delivered to the user during CPAP therapy.

Compatibility	ResMed Airsense™ S9 & 10
Classification	Class I according to (EU) 2017/745, Annex VIII
Dimensions & Tolerances	51.0 × 35.0 mm
Color	White
Filter Type	Disposable
Outer Material	Non-woven polypropylene
Filter Media	Electrostatic non woven cloth
Packaging	12 pcs per pack
Filter Efficiency	99% (NaCl, 0.3 µm @ 5.33 cm/s, ±2%)
Airflow Resistance / Pressure Drop	Tested according to ISO 29463
Sterility	Non-sterile
Storage Conditions	Temperature: +5 °C to +40 °C Relative humidity: < 80%, non-condensing
Transport Conditions	Temperature: -20 °C to +50 °C
Shelf Life	5 years from the date of manufacture when stored in a sealed bag and protected from direct sunlight