


Air Filter for Löwenstein prismaLINE

Item No.: HC-161



Product Function:

CPAP filters are designed to reduce airborne contaminants such as dust, pollen, and microorganisms from the intake air, ensuring that clean and filtered air is delivered to the user during CPAP therapy.

Compatibility	Löwenstein prismaLine
Classification	Class I according to (EU) 2017/745, Annex VIII
Dimensions & Tolerances	60.0 × 36.0 mm × 6.0 mm
Color	Black
Filter Type	Disposable
Outer Material	Foam
Filter Media	Polyurethane foam
Packaging	2 pcs per pack
Filter Efficiency	99.9% (NaCl, 0.3 µm @ 5.33 cm/s, ±2%)
Airflow Resistance / Pressure Drop	Tested according to ISO 29463
Sterility	Non-sterile
Storage Conditions	Temperature: +5 °C to +40 °C Relative humidity: < 80%, non-condensing
Transport Conditions	Temperature: -20 °C to +50 °C
Shelf Life	5 years from the date of manufacture when stored in a sealed bag and protected from direct sunlight