



## CNSAC Nasal/Oral Cannula for Pressure Sensor

**Highly-sensitive Nasal/Oral Cannula with excellent wearing comfort**

**For a reliable diagnosis of obstructive sleep apnea according to AASM**

CNSAC Nasal/Oral Cannulas guarantee a high wearing comfort and a stable Luer-Lock connection for measuring air flow and snoring during sleep screening and sleep studies (PSG). Soft nasal prongs provide maximum wearing comfort, especially for long wearing periods.

### Summary advantages:

- ✓ Innovative design
- ✓ Easy application
- ✓ Higher comfort for patients
- ✓ Universal Luer-lock connection

## Accessories and Consumables for Cardiology, Neurology and Sleep Medicine

Specification Data					
Name	Material	Color	Diameter	Length	Inner diameter
Prong support tube	Medical grade PVC	Transparent	3.35 mm	50 cm	/
Main tube	Medical grade PVC	Transparent	4.3 mm	60 cm	/
Tubing connector	Medical grade PVC	Transparent	/	/	1.9 mm
Adult Prong dipped, Nasal oral	Medical grade PVC	Transparent	42.86 x 35.35 mm	/	1.31 mm
Oxygen tube 3-chanel lumen	Medical grade PVC	Transparent	/	24.47 mm	/

Specification Data	
Product weight:	15.35 g
Sterile	NO
Operation ambient temperature	15 °C- 35 °C
Storage condition	-40 °C – 60 °C ≤80% humidity

### CNSAC MedShop GmbH

#### CONTACT

- Telephone +49 (931) 359 094-65  
(Mo.-Fr.: 08:00-19:00) +43 (274) 290 013-400
- Fax +49 (931) 359 094-988  
+43 (274) 290 013-410
- E-Mail [info@cnsac-medshop.com](mailto:info@cnsac-medshop.com)
- Online [www.cnsac-medshop.com](http://www.cnsac-medshop.com)

